



## **ONDA**

The Coolwaves™  
Revolution of Body Contouring.  
Smart and Enhanced Microwaves  
that Excel at All Points:

- Localized fat
- Cellulite
- Skin tightening

**DEKA**  
Innate Ability

Our smart  
microwaves  
are so cool,  
that your patients  
can stay cool  
too!



ONDA



## DEKA microwaves are not just smart. They are cool!

Here at DEKA, we have come up with yet another groundbreaking technological innovation: ONDA, the only system that uses Coolwaves™, special microwaves that selectively target fat cells to reduce localized fat deposits in a safe, effective and non-invasive way.

DEKA's R&D department has developed a smart handpiece (patent pending) that guides the electromagnetic field to the lipocytes in order to destroy them, while the integrated cooling system preserves the skin, ensuring patient safety and comfort.

“ In 1999 DEKA pioneered minimally invasive treatment of localized fat deposits by introducing Laserlipolysis® which achieved worldwide success. Today DEKA has taken one step further. DEKA's unceasing search for new sources and applications has led to the invention of a completely non-invasive system, which specifically destroys the lipocyte membrane, just as the laser did. In just a few sessions (usually 4-5), the Coolwaves™ in the ONDA system are extremely effective at reducing the localized accumulations of fat in such areas as the abdomen, back, trochanter and lower limbs. And that's not all! Coolwaves™ make it possible to treat cellulite even in the most advanced stages and skin laxity. This is a completely painless procedure, with no side effects, which is safe and comfortable for my patients, who were enthusiastic about this innovative treatment. ”

Paolo Bonan, M.D.  
ESLD Key Officer Education  
In Charge of Laser Cutaneous Cosmetic & Plastic Surgery Unit  
Villa Donatello Clinic, Florence - Italy



## The First Body Contouring System that Excels at All Points

DEKA's new microwave source was developed to provide a non-invasive system, in order to effectively intervene on the main blemishes requiring body remodeling:

- **Localized fat deposits:** The deep handpiece enables Coolwaves™ to penetrate deeply, disrupting the membranes of subcutaneous fat cells. The lysate is subsequently drained off by stimulating physiological metabolic processes (macrophage lipolysis). The cooling system integrated into the handpieces ensures that lipocytes are broken down in total safety for the skin and maximum patient comfort.
- **Cellulite:** The shallow handpiece effectively targets the connective septa between the adipose lobules in the more advanced stages of cellulite (visible in the typical "orange peel" effect).
- **Skin tightening:** The ONDA system's special microwaves induce immediate shrinkage of the collagen fibres in the dermis and stimulate the production of new collagen. This makes tissues more compact and toned.



### Onda's PLUS

Coolwaves™	ONDA is the only system that uses special microwaves to safely target subcutaneous fat with maximum patient comfort.
All-In-One Solution	A single, effective solution for localized fat deposits, cellulite and skin laxity.
Smart Handpieces	A system of LED lights enable smart handpieces (patent pending) to provide an intuitive guide to the operator during the procedure, which ensures that treatment is both safe and effective.
Skin Cooling	Integrated contact skin cooling system.
Simple & User Friendly	DEKA's R&D team gave ONDA an ergonomic, attractive design featuring a compact structure with modern lines; intuitive simple software that helps the operator perform the treatment, while the cooled, ergonomic handpieces ensure the treatment is carried out properly and provide the operator with warnings if any procedures are performed incorrectly.
Great Benefits	A rapid Return On Investment (ROI) and total patient satisfaction will increase your business!

### Specially designed smart handpieces: Deep & Shallow.



The internal mechanics of ONDA handpieces are specially designed to emit Coolwaves™ in a controlled fashion, focusing energy only where it is needed (for example in the subcutaneous adipose tissue in the event of localized fat deposits), preserving the adjacent tissues. The two handpieces differ in their depth of action, each designed for different body areas or types of blemish.



## Coolwaves™: a New Concept of Energy

Microwaves are part of the large family of electromagnetic waves [frequency range: 1 - 300 GHz]. Today they are common in many fields of everyday life including the medical field, for diagnostic or treatment purposes.

**However, not all microwaves are the same. No other device emits microwaves in the same way as the ONDA system's smart handpieces!**

Microwaves cause the molecules that absorb them to oscillate and vibrate, thus heating up the medium that contains the molecules.



# 80%

Coolwaves™ energy targeting the fat cells



# 20%



Secondary superficial thermal effect

The ONDA system handpieces generate microwaves at 2.45 GHz, a frequency that is preferentially absorbed by fat molecules rather than by water molecules. The conductivity of the epidermis and the dermis, which contain plenty of water but not of fat, is therefore much greater than that of the hypodermis. Therefore, when the Coolwaves™ penetrate the skin, they pass straight through the top layers without overheating them (using about 20% of the energy from the microwaves), while they concentrate their effective action on the subcutaneous fat that absorbs them more (about 80% of the energy). By contrast, traditional radio frequency systems (frequency range 0.1 - 40 MHz) act only on the surface due to their high absorption by the water molecules. This means that they not only fail to penetrate deep into the body to treat fat, but they also risk damaging the skin.




# Coolwaves™ & Smart Handpieces: Guarantee of Safety and Effectiveness

For further safety and improved patient comfort, the handpieces are equipped with an integrated cooling system that acts as a barrier between the Coolwaves™ delivery source and the skin. This creates a thermal cushion which protects the epidermis and the dermis, allowing the action to be concentrated deep down in the skin on the targeted fat. Coolwaves™ therefore do not risk causing hot spots on the skin, unlike the electrodes used in radiofrequency systems.




(A): ONDA smart handpieces (patent pending) are designed so that only 20% of the energy goes to the dermis, and this is in any case counterbalanced by the integrated cooling system that annuls the effects of such heat. The remaining 80% of the energy penetrates into the fat, acting effectively on the lipocytes. (B): the situation is quite different with RF handpieces. For Coolwaves, the conductivity of the outermost layers of skin is at least 3.5 times higher than that of the commonly-used RF irradiation systems in aesthetic medicine (Graph C). That means that most of the RF energy gets stuck in the epidermis and dermis, heating them up to such an extent that there is a risk of tissue damage. Moreover, as the RF energy remains close to the surface, it fails to reach the hypodermis where the fat cells are located, and whose membranes must be broken in order for the treatment to be effective.




(A) Thermographic image of ex-vivo tissue, seen in section, treated with an ONDA handpiece. A "hot zone" can clearly be seen just below the top layers of the skin, which stay cold and are therefore displayed in blue. The image next to it (B) corresponds to a similar ex-vivo tissue treated with bipolar radiofrequency. The tip effect where the electrodes and the skin come into contact clearly creates a strong rise in temperature on the surface, while deeper the tissue remains cold, which proves that it is not possible to reach the correct treatment temperatures in the subcutaneous fat.



## Histological & Clinical Study




A




B

H&E histological images of human abdominal tissue. (A): control. (B): Sample from the same patient immediately after treatment with ONDA Coolwaves™. When we compare the two images, it is clear that the epidermis in the treated sample has not undergone any alteration compared with the control. Moreover, in image (B), the collagen in the dermis looks more eosinophilous (pink): heat causes shrinking or tightening of collagen, resulting in greater eosinophilia, due to the higher collagen concentration, and is always accompanied by the presence of lighter-coloured spaces clearly visible in H&E, freed up by the aggregating collagen.



A



B

Histological images of tissue with human abdominal fat. (A): Control (magnification x10). (B): Sample from the same patient immediately after treatment with ONDA Coolwaves™ (magnification x20). Image (B) clearly shows the ruptured lipocytes and initial hyperaemia with dilatation of the blood vessels.

*Courtesy of Prof. R. Perrotta, M.D. and M.S. Tarico, M.D., Catania - Italy*

The Coolwaves™ technology is quick, easy to learn and can be put immediately into action by different categories of staff (doctors, medical assistants, nurses).




A



B

Before and after 1 session.

*Courtesy of Prof. R. Perrotta, M.D. and M.S. Tarico, M.D., Catania - Italy*



A



B

Before and after 1 session.

*Courtesy of Idaga Centre, Ondarroa - Spain*



## Simple, User-Friendly and Intuitive Software

**Built-in  
database protocol:  
expert  
from the start!**

### Graphic User Interface (GUI)

A new simple, intuitive touch-screen interface helps the operator select the most appropriate protocol, based on the specific features of each individual patient. The large LCD Touch Screen offers a quick and easy selection of the operating parameters.

From an operational viewpoint the procedure is very simple:

- The handpieces are automatically recognized by the system.
- The simplified software and integrated database permit rapid choice of parameters for any kind of treatment and body area.

- An accumulator displays the dose delivered during the treatment and an acoustic signal warns the user when the desired end-point has been reached. At this point the operator can treat another area.



multimedia



user-friendly



design



simplicity



## Technical Data

Onda	
Source	Microwaves - Coolwaves™
Frequency	2.45 GHz
Power	250 W
Smart Handpieces (Patent Pending)	Deep: Ø 6.6 cm - 35 cm <sup>2</sup> Shallow: Ø 5.6 cm - 25 cm <sup>2</sup> LED system control to provide an intuitive guide to the operator during the procedure
Handpiece Cooling	Integrated Skin Cooler
Emission Control	Fingerswitch
Graphical User Interface	10.4" Colour Display Touch Screen 10.4"
Data Base	Integrate tutorial with treatment protocol
Dimensions	394 (W) mm x 1075 (H) mm x 738 (D) mm
Weight	< 60 kg
Electrical Requirements	100-230 Vac, 50-60 Hz, 1500 VA



This brochure is not intended for the market of USA.

**DEKA**  
Innate Ability

Follow us on



[www.dekalaser.com](http://www.dekalaser.com)

**CE**  
0123 PENDING

onda

**E.E.**  
-El.En. Group

Dealer stamp

**DEKA M.E.L.A. s.r.l.**

Via Baldanzese,17  
50041 Calenzano (FI) - Italy  
Tel. +39 055 8874942 - Fax +39 055 8832884

**DEKA Innate Ability**

A spin-off of the El.En. Group, DEKA is a world-class leader in the design and manufacture of lasers and light sources for applications in the medical field. DEKA markets its devices in more than 80 countries throughout an extensive network of international distributors as well as direct offices in Italy, France, Japan and USA. Excellence is the hallmark of DEKA's experience and recognition garnered in the sphere of R&D in over thirty years of activity. Quality, innovation and technological excellence place DEKA and its products in a unique and distinguished position in the global arena. DEKA manufactures laser devices in compliance with the specifications of Directive 93/42/EEC and its quality assurance system is in accordance with the ISO 9001 and ISO 13485 standards.

DEKA M.E.L.A. s.r.l. - All rights reserved - In order to improve its products the company reserves the right to modify these specifications without prior notice. Document Reserved for Health Professionals Only.